



Disposable Suction Tubes for Ear, Nose, Troath & Sinus



Product Overview

Suction Tube for Ear Disposable

Suction tubes for Ear are designed for suction of fluid and other media from the ear. The malleable shaft allows for additional bending. Fine suction devices for even the narrowest approaches.

Design Benefits:

- Very stable suction tube for retraction and suction.
- The curved angle is specially designed for microscopic examination and surgery.
- Ideal for precise dissections.
- Quick connector Suction tip handle is easy to change during surgery.
- Reduces unnecessary infection risks.



Reference	Description	Unit/Box
FI153080	Suction Tube for Ear-Disposable (3.0Dia.x80.0)mm	50 pcs
FI152580	Suction Tube for Ear-Disposable (2.5Dia.x80.0)mm	50 pcs
FI1525100	Suction Tube for Ear-Disposable (2.5Dia.x100.0)mm	50 pcs
FI152060	Suction Tube for Ear-Disposable (2.0Dia.x60.0)mm	50 pcs
FI152070	Suction Tube for Ear-Disposable (2.0Dia.x70.0)mm	50 pcs
FI152080	Suction Tube for Ear-Disposable (2.0Dia.x80.0)mm	50 pcs
FI10110	Suction Tube for Ear-Disposable (2.0Dia.x100.0)mm	50 pcs
FI151670	Suction Tube for Ear-Disposable (1.6Dia.x70.0)mm	50 pcs
FI10111	Suction Tube for Ear-Disposable (1.6Dia.x100.0)mm	50 pcs
FI10112	Suction Tube for Ear-Disposable (1.4Dia.x100.0)mm	50 pcs
FI151480	Suction Tube for Ear-Disposable (1.4Dia.x80.0)mm	50 pcs
FI151470	Suction Tube for Ear-Disposable (1.4Dia.x70.0)mm	50 pcs
FI151460	Suction Tube for Ear-Disposable (1.4Dia.x60.0)mm	50 pcs
FI10113	Suction Tube for Ear-Disposable(1.2Dia.x100.0)mm	50 pcs
FI151270	Suction Tube for Ear-Disposable (1.2Dia.x70.0)mm	50 pcs
FI10114	Suction Tube for Ear-Disposable (1.0Dia.x100.0)mm	50 pcs
FI151080	Suction Tube for Ear-Disposable (1.0Dia.x80.0)mm	50 pcs
FI151070	Suction Tube for Ear-Disposable (1.0Dia.x70.0)mm	50 pcs
FI151060	Suction Tube for Ear-Disposable (1.0Dia.x60.0)mm	50 pcs
FI10115	Suction Tube for Ear-Disposable (0.8Dia.x100.0)mm	50 pcs
FI150870	Suction Tube for Ear-Disposable (0.8Dia.x70.0)mm	50 pcs
FI10116	Suction Tube for Ear-Disposable (0.7Dia.x100.0)mm	50 pcs
FI150780	Suction Tube for Ear-Disposable (0.7Dia.x80.0)mm	50 pcs
FI150770	Suction Tube for Ear-Disposable (0.7Dia.x70.0)mm	50 pcs
FI150760	Suction Tube for Ear-Disposable (0.7Dia.x60.0)mm	50 pcs
FI150670	Suction Tube for Ear-Disposable (0.6Dia.x70.0)mm	50 pcs
FI10117	Suction Tube for Ear-Disposable (0.4Dia.x100.0)mm	50 pcs
FI150480	Suction Tube for Ear-Disposable (0.4Dia.x80.0)mm	50 pcs
FI150470	Suction Tube for Ear-Disposable (0.4Dia.x70.0)mm	50 pcs
FI150460	Suction Tube for Ear-Disposable (0.4Dia.x60.0)mm	50 pcs
FI10120	Suction Tube for Ear-Disposable-SET*	8 pcs

* F10110, F10111, F10112, F10113, F10114, F10115, F10116, F10117

Suction Tube for Nose Disposable

Suction tubes for Nose are used during Rhinology and evacuates liquid from nasal cavities.

Design Benefits:

- Very stable suction tube for retraction and suction.
- The curved angle is specially designed for microscopic examination and surgery.
- Ideal for precise dissections.
- Quick connector Suction tip Handle is easy to change during surgery.
- Reduces unnecessary infection risks.



Reference	Description	Unit/Box
FI1220100	Suction Tube for Nose-Disposable (2.0Dia.x100.0)mm	25 pcs
FI1220120	Suction Tube for Nose-Disposable (2.0Dia.x120.0)mm	25 pcs
FI1220140	Suction Tube for Nose-Disposable (2.0Dia.x140.0)mm	25 pcs
FI1220150	Suction Tube for Nose-Disposable (2.0Dia.x150.0)mm	25 pcs
FI1220160	Suction Tube for Nose-Disposable (2.0Dia.x160.0)mm	25 pcs
FI1225120	Suction Tube for Nose-Disposable (2.5Dia.x120.0)mm	25 pcs
FI1230100	Suction Tube for Nose-Disposable (3.0Dia.x100.0)mm	25 pcs
FI1230120	Suction Tube for Nose-Disposable (3.0Dia.x120.0)mm	25 pcs
FI1230140	Suction Tube for Nose-Disposable (3.0Dia.x140.0)mm	25 pcs
FI1230150	Suction Tube for Nose-Disposable (3.0Dia.x150.0)mm	25 pcs
FI1230160	Suction Tube for Nose-Disposable (3.0Dia.x160.0)mm	25 pcs
FI1235120	Suction Tube for Nose-Disposable (3.5Dia.x120.0)mm	25 pcs
FI1235140	Suction Tube for Nose-Disposable (3.5Dia.x140.0)mm	25 pcs
FI1235160	Suction Tube for Nose-Disposable (3.5Dia.x160.0)mm	25 pcs
FI1240100	Suction Tube for Nose-Disposable (4.0Dia.x100.0)mm	25 pcs
FI1240120	Suction Tube for Nose-Disposable (4.0Dia.x120.0)mm	25 pcs
FI1240140	Suction Tube for Nose-Disposable (4.0Dia.x140.0)mm	25 pcs
FI1240150	Suction Tube for Nose-Disposable (4.0Dia.x150.0)mm	25 pcs
FI1240160	Suction Tube for Nose-Disposable (4.0Dia.x160.0)mm	25 pcs
FI1250120	Suction Tube for Nose-Disposable (5.0Dia.x120.0)mm	25 pcs
FI1250140	Suction Tube for Nose-Disposable (5.0Dia.x140.0)mm	25 pcs
FI1250160	Suction Tube for Nose-Disposable (5.0Dia.x160.0)mm	25 pcs

Suction Tube for Throat & Sinus Disposable

Suction tubes for Throat and Sinus are used to remove excess amount of foreign body material from throat and sinus.

Design Benefits:

- Very stable suction tube for retraction and suction.
- The shaft is slightly curved
- Reduces unnecessary infection risks.

Reference	Description	Unit/Box
FI1620250	Suction Tube for Throat-Disposable (2.0Dia.x250.0)mm	50 pcs
FI1630250	Suction Tube for Throat-Disposable (3.0Dia.x250.0)mm	50 pcs
FI1825150	Suction Tube for Sinus-Disposable (2.5Dia.x150.0)mm	50 pcs
FI1830150	Suction Tube for Sinus-Disposable (3.0Dia.x150.0)mm	50 pcs
FI1840150	Suction Tube for Sinus-Disposable (4.0Dia.x150.0)mm	50 pcs



Suction Tube for Throat and its tip close-up
(Material: SS and Polyethylene)



Suction Tube for Sinus and its tip close-up
(Material: SS and Polyethylene)

REF	Description	Unit/Box
EI0968	Suction Controller, Lenght : 39mm	50 pcs



During the suction procedure it is define in which level suction will perform for that precise suction level suction controller is require. Controls the level of suction via the thumb vent mechanism.

Design Benefits:

- Assemble with the desired gauge of suction tube for a precise and controllable procedure / Easy to handle

Material: Polyethylene