

# ELECTROSURGICAL ACCESSORIES

TECHNICAL INFORMATION



# INDEX

---

- Page 06: *Power Supply Cable*
- Page 09: *Surgical pencil*
- Page 12: *Foot Switch*
- Page 16: *Neutral*
- Page 23: *Connection*
- Page 29: *Electrode*
- Page 32: *Adapters*
- Page 34: *Bipolar/Monopolar*
- Page 41: *Evacuation Kit*
- Page 45: *Units support*

# CAPTIONS

## Compatibility:

- + = optional
- x = not compatible
- # = to verify
- N = Quantity Supplied



Packaging:  
quantity of articles contained



Autoclavable at 134°



Not reusable



Sterilized with ethylene

IP67

The item is completely dust-proof  
and can withstand a 30-minute  
immersion in standing water




Anti explosion



100 Sterilization cycles



30 Sterilization cycles



A black and white photograph showing a close-up of a person's hands wearing white nitrile gloves. The person is holding a clear plastic container with a lid, which appears to be a medical or laboratory sample. They are also wearing a white surgical mask. The background is slightly blurred.

AC  
CÈS  
SÖES  
RIES

A black and white photograph of two medical professionals, likely surgeons, wearing surgical caps and masks. They are positioned in an operating room, with one surgeon leaning forward over a patient and the other standing beside them. A large, circular surgical light fixture is visible above them. The scene conveys a sense of medical precision and teamwork.

**POWER  
SUPPLY  
CABLE**

<b>Description</b>	Power supply cable SIE-IEC
<b>Cable Length</b>	5 m (196,85 in)

00100.01



<b>Description</b>	Power supply cable SIE-IEC
<b>Cable Length</b>	2 m (78,74 in)

00100.03




<b>Description</b>	Power supply cable IT-IEC
<b>Cable Length</b>	2 m (78,74 in)

00100.00



**Description** Power supply cable USA-IEC  
**Cable Length** 2 m (78.74 in)

00100.04



5 m



2 m




2 m



2 m

00100.05




2 m

<b>Description</b>	Power supply cable UK-IEC
<b>Cable Length</b>	2 m (78,74 in)



00100.07



2 m

<b>Description</b>	Power supply cable BR-IEC
<b>Cable Length</b>	2 m (78,74 in)



00100.09



2 m

<b>Description</b>	Power supply cable AU/CN-IEC
<b>Cable Length</b>	2 m (78,74 in)



pcs: 1

00100.10



2 m


<b>Description</b>	Power supply cable JP-IEC
<b>Cable Length</b>	2 m (78,74 in)



A black and white photograph showing two hands wearing white medical gloves. The hands are positioned as if they are about to use a surgical pencil, with fingers gripping the handle and the lead tip pointing downwards. The background is a solid dark color.

# SURGICAL PENCIL

## 00201.02



### Description

Handle for microsurgical needle autoclavable

### Cable Length

2 m (78,74 in)

### Size

10x1x1 cm (3,93x0,39x0,39 in)

### Weight

0,05 kg (0,11lb)

### Electrode Compatibility

ø0,8 mm (0,031 in)

### Sterilization

Autoclavable

### Sterilization cycles

100

### Compatibility

500500.L



pcs: 1

## 00205.00



### Description

Handle with switches

3 m (118,11 in)

### Size

18,5x1,6x1,7 cm (7,2x0,62x0,66 in)

### Weight

0,03kg (0,06 lb)

### Electrode Compatibility

ø2,4 mm (0,94 in)

### Sterilization

Autoclavable

### Sterilization cycles

100



pcs: 1

## 00205.40



### Description

Handle with switches

3 m (118,11 in)

### Size

18,5x1,6x1,7 cm (7,2x0,62x0,66 in)

### Weight

0,03kg (0,06 lb)

### Electrode Compatibility

ø4 mm (1,15 in)

### Sterilization

Autoclavable

### Sterilization cycles

100

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x x	+ +	+ +	+ +	+ +	+ +	1 1	1 1	1 1	1 1	1 1	1 1	1 1	1 1	x x	x x



pcs: 1



# 00206.00

<b>Description</b>	Monopolar Handpiece
<b>Cable Length</b>	3 m (118,11 in)
<b>Size</b>	18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
<b>Weight</b>	0,03kg (0,06 lb)
<b>Electrode Compatibility</b>	ø2,4 mm (0,94 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100



pcs: 1

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
1 1	+	+	+	+	+	+	+	+	+	+	+	+	x x	x x

<b>Description</b>	Monopolar Handpiece
<b>Cable Length</b>	3m (118,11 in)
<b>Size</b>	18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
<b>Weight</b>	0,03 kg (0,06 lb)
<b>Electrode Compatibility</b>	ø4 mm (1,15 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100



pcs: 1

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+	+	+	+	+	+	+	+	+	+	+	+	+	x x	x x

# 00206.40 S



pcs: 1

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+	+	+	+	+	+	+	+	+	+	+	+	+	x x	x x

**Description**


<b>Description</b>	Handpiece for argon electrodes reusable
<b>Cable Length</b>	3,5 m (137,8 in)
<b>Size</b>	18,5x1,6x1,7 cm (7,2x0,62x0,66 in)
<b>Weight</b>	0,03kg (0,06lb)
<b>Electrode Compatibility</b>	ø4 mm (1,15 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100



pcs: 1

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x x	x x	x x	x x	x x	x x	+	+	+	+	+	+	+	x x	x x

# 00420.00

**LUER LOCK**

3xMØ4 mm

**Description**

<b>Description</b>	Disposable handle with finger switches (F4797)
<b>Size</b>	14x2x1 cm (5,5x0,78x0,39 in)
<b>Cable Length</b>	3m (118,11in)
<b>Weight</b>	0,07 kg (0,15 lb)
<b>Electrode Compatibility</b>	ø2,4 mm (0,094 in)
<b>Sterilization</b>	Sterile EO
<b>Sterilization cycles</b>	0



pcs: 1

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x x	+	+	+	+	+	+	+	+	+	5 5	5	5	x x	x x

# 755VL



STERILE EO



# FOOT SWITCH

00304.00

**Description**  
**Cable Length**  
**Size**  
**Weight**

Water-proof foot switch  
2 m (78,75 in)  
11x3,5x9 cm (7,2x0,62x0,66 in)  
0,3 kg (0,66 lb)

50 D	1	80 D	1	80	1	120	1	160	1	200	1
TOUCH 200	1	FLASH 120	1	FLASH 160	1	FLASH 200	1	x 300 HP	x x x x	400 HP	1

TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
# # # x x x x + + + +

300 HP  
400 HP  
TOUCH HP  
TOUCH HP (Argon)  
SURTRON EVAC  
SURTRON ABC



pcs: 1

**Description**  
**Cable Length**  
**Size**  
**Weight**

Single water-proof foot switch  
(3 pins)  
5 m (196,85 in)  
11x3,5x9 cm (7,2x0,62x0,66 in)  
0,3 kg (0,66 lb)

50 D	x	80 D	x	80	x	120	x	160	x	200	x
TOUCH 200		FLASH 120		FLASH 160		FLASH 200		x 300 HP	x x x x	400 HP	

TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
# # # x x x x + + + +

300 HP  
400 HP  
TOUCH HP  
TOUCH HP (Argon)  
SURTRON EVAC  
SURTRON ABC



pcs: 1

**Description**  
**Cable Length**  
**Size**  
**Weight**

Water-proof foot switch  
3 m (118,11 in)  
11x3,5x9 cm (7,2x0,62x0,66 in)  
0,3 kg (0,66 lb)

+ 50 D	+	+ 80 D	+	+ 80	+	+ 120	+	+ 160	+	+ 200	+
TOUCH 200		FLASH 120		FLASH 160		FLASH 200		x 300 HP	x x x x	400 HP	

TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
x x x x + + + +

300 HP  
400 HP  
TOUCH HP  
TOUCH HP (Argon)  
SURTRON EVAC  
SURTRON ABC



pcs: 1



WATER PROOF

IP67

00304.01

**Description**  
**Cable Length**  
**Size**  
**Weight**

Single water-proof foot switch  
(3 pins)  
5 m (196,85 in)  
11x3,5x9 cm (7,2x0,62x0,66 in)  
0,3 kg (0,66 lb)

50 D	x	80 D	x	80	x	120	x	160	x	200	x
TOUCH 200		FLASH 120		FLASH 160		FLASH 200		x 300 HP	x x x x	400 HP	

TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
# # # x x x x + + + +

300 HP  
400 HP  
TOUCH HP  
TOUCH HP (Argon)  
SURTRON EVAC  
SURTRON ABC



pcs: 1



WATER PROOF

IP67

00304.04

**Description**  
**Cable Length**  
**Size**  
**Weight**

Water-proof foot switch  
3 m (118,11 in)  
11x3,5x9 cm (7,2x0,62x0,66 in)  
0,3 kg (0,66 lb)

+ 50 D	+	+ 80 D	+	+ 80	+	+ 120	+	+ 160	+	+ 200	+
TOUCH 200		FLASH 120		FLASH 160		FLASH 200		x 300 HP	x x x x	400 HP	

TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
x x x x + + + +

300 HP  
400 HP  
TOUCH HP  
TOUCH HP (Argon)  
SURTRON EVAC  
SURTRON ABC




pcs: 1



WATER PROOF

IP67

## 00304.PO



**Description** Protection for single foot switch  
**Cable Length** NO  
**Size** 12x10x14 cm (4,7x3,9x1,64 in)  
**Weight** 0,65 kg (0,143 lb)

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+	+	+	+	+	+	+	+	+	+	x	x	x	x	#	#



## 00305.03



WATER PROOF IP67

**Description** Double water proof foot switch  
**Cable Length** 5 m (196,85 in)  
**Size** 35x6x18 cm (13,78x2,36x7 in)  
**Weight** 1,3 kg (2,86lb)

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x	x	+	+	+	+	+	+	+	+	x	x	x	x	x	x



## 00305.04



WATER PROOF IP67

**Description** Waterproof double foot switch  
**Cable Length** 5 m (196,85 in)  
**Size** 35x6x18 cm (13,78x2,36x7 in)  
**Weight** 1,3 kg (2,86lb)

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x	x	x	x	x	x	x	x	x	x	2	2	x	x	x	x



*00305.06*

<b>Description</b>	Waterproof double foot switch
<b>Cable Length</b>	5 m (196,85 in)
<b>Size</b>	35x6x18 cm (13,78x2,36x7 in)
<b>Weight</b>	1.3 kg (2,86lb)

50 D	TOUCH 200
x x 80 D	FLASH 120
80	FLASH 160
x x 120	FLASH 200
160	300 HP
x x 200	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



<b>Description</b>	Waterproof double foot switch with selector
<b>Cable Length</b>	5 m (196,85 in)
<b>Size</b>	35x6x18 cm (13,78x2,36x7 in)
<b>Weight</b>	1.3 kg (2,86lb)

50 D	TOUCH 200
x x 80 D	FLASH 120
80	FLASH 160
x x 120	FLASH 200
160	300 HP
x x 200	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



<b>Description</b>	Waterproof double foot switch with selector
<b>Cable Length</b>	5 m (196,85 in)
<b>Size</b>	35x6x18 cm (13,78x2,36x7 in)
<b>Weight</b>	1.3 kg (2,86lb)

50 D	TOUCH 200
x x 80 D	FLASH 120
80	FLASH 160
x x 120	FLASH 200
160	300 HP
x x 200	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



WATER PROOF IP67

*00306.03*

<b>Description</b>	Waterproof double foot switch with selector
<b>Cable Length</b>	5 m (196,85 in)
<b>Size</b>	35x6x18 cm (13,78x2,36x7 in)
<b>Weight</b>	1.3 kg (2,86lb)

50 D	TOUCH 200
x x 80 D	FLASH 120
80	FLASH 160
x x 120	FLASH 200
160	300 HP
x x 200	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



WATER PROOF IP67

*00306.06*

<b>Description</b>	Waterproof double foot switch with selector
<b>Cable Length</b>	5 m (196,85 in)
<b>Size</b>	35x6x18 cm (13,78x2,36x7 in)
<b>Weight</b>	1.3 kg (2,86lb)

50 D	TOUCH 200
x x 80 D	FLASH 120
80	FLASH 160
x x 120	FLASH 200
160	300 HP
x x 200	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC




WATER PROOF IP67



NEU  
TRAL

<b>Description</b>	Rod neutral electrode with cable
<b>Cable Length</b>	1,8 m (70,86 in)
<b>Size</b>	12xØ1,8 cm (4,72x0,7 in)
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO

00400.00



+ 50 D	+ 80 D	+ 80
+ + 80	+ + 120	
+ + 120	+ + 160	
+ + 160	+ + 200	
+ + 200		
	+ TOUCH 200	
	+ FLASH 120	
	+ FLASH 160	
	+ FLASH 200	
	+ 300 HP	
	+ 400 HP	
	+ TOUCH HP	
	+ TOUCH HP (Argon)	
	+ SURTRON EVAC	
	+ SURTRON ABC	




pcs: 1

MAX. 120W

Integrated connection

<b>Description</b>	Steel neutral electrode with cable
<b>Cable Length</b>	3 m (118,11 in)
<b>Size</b>	16x12 cm (6,3x4,7 in)
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO

00401.00



+ 50 D	+ 80 D	+ 80
+ + 80	+ + 120	
+ + 120	+ + 160	
+ + 160	+ + 200	
+ + 200		
	+ TOUCH 200	
	+ FLASH 120	
	+ FLASH 160	
	+ FLASH 200	
	+ 300 HP	
	+ 400 HP	
	+ TOUCH HP	
	+ TOUCH HP (Argon)	
	+ SURTRON EVAC	
	+ SURTRON ABC	




pcs: 1

Integrated connection

<b>Description</b>	Steel neutral electrode with cable
<b>Cable Length</b>	3 m (118,11 in)
<b>Size</b>	24x16 cm (9,4x6,3 in)
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO

00401.01




+ 50 D	+ 80 D	+ 80
+ + 80	+ + 120	
+ + 120	+ + 160	
+ + 160	+ + 200	
+ + 200		
	+ TOUCH 200	
	+ FLASH 120	
	+ FLASH 160	
	+ FLASH 200	
	+ 300 HP	
	+ 400 HP	
	+ TOUCH HP	
	+ TOUCH HP (Argon)	
	+ SURTRON EVAC	
	+ SURTRON ABC	



pcs: 1

Integrated connection

## 00401.02




<b>Description</b>	Steel neutral electrode autoclavable
<b>Size WxHxD</b>	16x12x0,4 cm (6,3x4,7x0,15 in)
<b>Cable Length</b>	3 m (118,11 in)
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	
+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	SURTRON ABC



## 00401.03



<b>Description</b>	Steel neutral electrode autoclavable
<b>Size WxHxD</b>	24x16x0,4 cm (9,4x6,3x0,15 in)
<b>Cable Length</b>	3 m (118,11 in)
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	
+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	SURTRON ABC



## 00401.10



<b>Description</b>	Steel neutral electrode FLEX
<b>Size WxHxD</b>	12x21x0,2 cm (4,7x8,2x0,078 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,01 Kg (0,02lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100
<b>Compatibility</b>	00404.08, 00404.09, 00404.10, 00404.11


50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)	SURTRON EVAC	
+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	+ + + + + + + + + + + + + +	SURTRON ABC



# 00401.11

<b>Description</b>	Steel neutral electrode FLEX with cable
<b>Size WxHxD</b>	12x21x0,2 cm (4,7x8,2x0,078 in)
<b>Cable Length</b>	3m (118,11 in)
<b>Weight</b>	0,01 Kg (0,02lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO

50 D  
+ 80 D  
+ + 80  
+ + + 120  
+ + + + 160  
+ + + + + 200  
+ + + + + + TOUCH 200  
+ + + + + + FLASH 120  
+ + + + + + FLASH 160  
+ + + + + + FLASH 200  
+ + + + + + 300 HP  
+ + + + + + 400 HP  
+ + + + + + TOUCH HP  
+ + + + + + TOUCH HP (Argon)  
+ + + + + + SURTRON EVAC  
+ + + + + + SURTRON ABC



Integrated connection

<b>Description</b>	Steel neutral electrode FLEX with cable autoclavable
<b>Size WxHxD</b>	12x21x0,2 cm (4,7x8,2x0,078 in)
<b>Cable Length</b>	3m (118,11 in)
<b>Weight</b>	0,01 Kg (0,02lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

50 D  
+ 80 D  
+ + 80  
+ + + 120  
+ + + + 160  
+ + + + + 200  
+ + + + + + TOUCH 200  
+ + + + + + FLASH 120  
+ + + + + + FLASH 160  
+ + + + + + FLASH 200  
+ + + + + + 300 HP  
+ + + + + + 400 HP  
+ + + + + + TOUCH HP  
+ + + + + + TOUCH HP (Argon)  
+ + + + + + SURTRON EVAC  
+ + + + + + SURTRON ABC



Integrated connection

134°C  
111°C  
134°C

CICLI 100 CYCLES

<b>Description</b>	Steel neutral electrode FLEX S
<b>Size WxHxD</b>	12x21x0,2 cm (4,7x8,2x0,078 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,01 Kg (0,02 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100
<b>Compatibility</b>	00404.08, 00404.09, 00404.10, 00404.11

50 D  
# + 80 D  
+ + 80  
+ + + 120  
+ + + + 160  
+ + + + + 200  
+ + + + + + TOUCH 200  
+ + + + + + FLASH 120  
+ + + + + + FLASH 160  
+ + + + + + FLASH 200  
+ + + + + + 300 HP  
+ + + + + + 400 HP  
+ + + + + + TOUCH HP  
+ + + + + + TOUCH HP (Argon)  
+ + + + + + SURTRON EVAC  
+ + + + + + SURTRON ABC



134°C  
111°C  
134°C

CICLI 100 CYCLES

## 00401.21



Integrated connection

**Description**

Steel neutral electrode FLEX S with cable

**Size WxHxD**

12x21x0,2 cm (4,7x8,2x0,078 in)

**Cable Length**

3 m (118,11 in)

**Weight**

0,01 Kg (0,02 lb)

**Sterilization**

NO

**Sterilization cycles**

NO

#	50 D	80 D	80	120	160	200	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
#	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x x	x x



pcs: 1

## 00401.22



Integrated connection



**Description**

Steel neutral electrode FLEX S with cable autoclavable

**Size WxHxD**

12x21x0,2 cm (4,7x8,2x0,078 in)

**Cable Length**

3 m (118,11 in)

**Weight**

0,01 kg (0,02 lb)

**Sterilization**

Autoclavable

**Sterilization cycles**

100

#	50 D	80 D	80	120	160	200	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
#	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	x x	x x



pcs: 1

## 00403.01



**Description**

Steel neutral electrode autoclavable

**Size WxHxD**

10x1,2x1,2 cm (3,934x0,47x0,47 in)

**Cable Length**

NO

**Weight**

0,1 kg (0,22 lb)

**Sterilization**

Autoclavable

**Sterilization cycles**

100

**Compatibility**

00404.02, 00404.04

1	x	x	80	80	120	160	200	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
1	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x x	x x	



pcs: 1

00403.02


<b>Description</b>	Rod neutral electrode
<b>Size WxHxD</b>	14x2x2 cm (5,5x0,78x0,78 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,4 Kg (0,88 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100
<b>Compatibility</b>	00404.02, 00404.04

50 D	80 D
+ x	x x
111	111
200	200
TOUCH 200	FLASH 200
FLASH 120	FLASH 120
FLASH 160	FLASH 160
FLASH 200	FLASH 200
300 HP	300 HP
400 HP	400 HP
TOUCH HP	TOUCH HP
TOUCH HP (Argon)	TOUCH HP (Argon)
SURTRON EVAC	SURTRON ABC
SURTRON ABC	SURTRON ABC




<b>Description</b>	Steel neutral electrode
<b>Size WxHxD</b>	16x12x0,02 cm (6,3x4,7x0,008 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100
<b>Compatibility</b>	00404.08, 00404.09, 00404.10, 00404.11

50 D	80 D
+ 1	1 1
111	111
200	200
TOUCH 200	FLASH 200
FLASH 120	FLASH 120
FLASH 160	FLASH 160
FLASH 200	FLASH 200
300 HP	300 HP
400 HP	400 HP
TOUCH HP	TOUCH HP
TOUCH HP (Argon)	TOUCH HP (Argon)
SURTRON EVAC	SURTRON ABC
SURTRON ABC	SURTRON ABC




<b>Description</b>	Steel neutral electrode
<b>Size WxHxD</b>	24x16x0,02 cm (9,4x6,3x0,008 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,1 kg (0,22 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100
<b>Compatibility</b>	00404.08, 00404.09, 00404.10, 00404.11

50 D	80 D
+ 1	1 1
111	111
200	200
TOUCH 200	FLASH 200
FLASH 120	FLASH 120
FLASH 160	FLASH 160
FLASH 200	FLASH 200
300 HP	300 HP
400 HP	400 HP
TOUCH HP	TOUCH HP
TOUCH HP (Argon)	TOUCH HP (Argon)
SURTRON EVAC	SURTRON ABC
SURTRON ABC	SURTRON ABC




F7920



2

F7930



<b>Description</b>	Disposable split neutral electrode (F7820)
<b>Size WxHxD</b>	20x10,5x0,2 cm (7,8x4,1x0,078 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,01 Kg (0,02 lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	00404.08, 00404.09, 00404.10, 00404.11



<b>Description</b>	Conductive rubber split neutral electrode
<b>Size WxHxD</b>	20x10,5x0,4 cm (7,8x4,1x0,15 in)
<b>Cable Length</b>	NO
<b>Weight</b>	0,1 Kg (0,2lb)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	00404.07




pcs: 1

# CON NECT ION



## 00402.00



<b>Description</b>	Monopolar Cable M4-F2.8
<b>Weight</b>	0,06 kg (0,13lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

134°C	CICLI 100 CYCLES	# 50 D	# 80 D	# 80	+ 120	+ 160	+ 200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	+ 300 HP	+ 400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
		# # 50 D	# # 80 D	# # 80	+ + 120	+ + 160	+ + 200	+ + TOUCH 200	+ + FLASH 120	+ + FLASH 160	+ + FLASH 200	+ + 300 HP	+ + 400 HP	+ + TOUCH HP	+ + TOUCH HP (Argon)	+ + SURTRON EVAC	+ + SURTRON ABC

pcs: 1

## 00402.01



<b>Description</b>	Monopolar Cable M4-F4
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

134°C	CICLI 100 CYCLES	# 50 D	# 80 D	# 80	+ 120	+ 160	+ 200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	+ 300 HP	+ 400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
		# # 50 D	# # 80 D	# # 80	+ + 120	+ + 160	+ + 200	+ + TOUCH 200	+ + FLASH 120	+ + FLASH 160	+ + FLASH 200	+ + 300 HP	+ + 400 HP	+ + TOUCH HP	+ + TOUCH HP (Argon)	+ + SURTRON EVAC	+ + SURTRON ABC

pcs: 1

## 00402.02




<b>Description</b>	Monopolar Cable M4-MP4
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

134°C	CICLI 100 CYCLES	# 50 D	# 80 D	# 80	+ 120	+ 160	+ 200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	+ 300 HP	+ 400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
		# # 50 D	# # 80 D	# # 80	+ + 120	+ + 160	+ + 200	+ + TOUCH 200	+ + FLASH 120	+ + FLASH 160	+ + FLASH 200	+ + 300 HP	+ + 400 HP	+ + TOUCH HP	+ + TOUCH HP (Argon)	+ + SURTRON EVAC	+ + SURTRON ABC

pcs: 1

## 00402.03



<b>Description</b>	Monopolar Cable M4-EU
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	100

134°C	CICLI 100 CYCLES	# 50 D	# 80 D	# 80	+ 120	+ 160	+ 200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	+ 300 HP	+ 400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
		# # 50 D	# # 80 D	# # 80	+ + 120	+ + 160	+ + 200	+ + TOUCH 200	+ + FLASH 120	+ + FLASH 160	+ + FLASH 200	+ + 300 HP	+ + 400 HP	+ + TOUCH HP	+ + TOUCH HP (Argon)	+ + SURTRON EVAC	+ + SURTRON ABC


pcs: 1

**Description** Monopolar Cable M4-F2.M  
**Weight** 0,06 kg (0,13 lb)  
**Cable Length** 3 m (118,11 in)  
**Sterilization** Autoclavable  
**Sterilization cycles** 100

50 D	80 D
# #	80
80	
+ +	120
120	
+ +	160
160	
+ +	200
200	
+	TOUCH 200
+	FLASH 120
+	FLASH 160
+	FLASH 200
+	300 HP
+	400 HP
+	TOUCH HP
+	TOUCH HP (Argon)
+	SURTRON EVAC
x x x x x x	SURTRON ABC



00402.04




**Description** Monopolar Cable M4-F2.M  
**Weight** 0,06 kg (0,13 lb)  
**Cable Length** 1,8 m (70,86 in)  
**Sterilization** NO  
**Sterilization cycles** NO

50 D	80 D
1 x	80
80	
+ +	120
120	
+ +	160
160	
+ +	200
200	
+	TOUCH 200
+	FLASH 120
+	FLASH 160
+	FLASH 200
+	300 HP
+	400 HP
+	TOUCH HP
+	TOUCH HP (Argon)
+	SURTRON EVAC
x x x x x x	SURTRON ABC



00404.02




**Description** Cable for rod neutral electrode  
**Autoclavable**  
**Cable Length** 3 m (118 in)  
**Unit Connection** 4 mm (0,157 in)  
**Sterilization** Autoclavable  
**Sterilization cycles** 100

50 D	80 D
+ x	80
80	
+ +	120
120	
+ +	160
160	
+ +	200
200	
+	TOUCH 200
+	FLASH 120
+	FLASH 160
+	FLASH 200
+	300 HP
+	400 HP
+	TOUCH HP
+	TOUCH HP (Argon)
+	SURTRON EVAC
x x x x x x	SURTRON ABC



00404.04




**Description** Cable for neutral plate F7930  
**Cable Length** 5m (196,85 in)  
**Unit Connection** 4 mm (0,157 in)  
**Accessory Connection** 2,8 mm (0,110 in)  
**Compatibility** F7930


# 50 D	80 D
+ +	80
+ +	120
+ +	160
+ +	200
+	TOUCH 200
+	FLASH 120
+	FLASH 160
+	FLASH 200
+	300 HP
+	400 HP
+	TOUCH HP
+	TOUCH HP (Argon)
+	SURTRON EVAC
x x x x x x	SURTRON ABC



00404.07



## 00404.08


**Description**

Cable for neutral electrode disposable type 5365/6429/FLEX

Weight 0,06 kg (0,13 lb)

Cable Length 3 m (118,11 in)

Sterilization NO


Sterilization cycles NO

Compatibility 0350, F790, 00401.10, 00401.20, 5365A, 6429A



pcs: 1

## 00404.09


**Description**

Cable for connection neutral disposable 5365-6429/FLEX

FLEX S AUTOCLAVE

Weight 0,06 kg (0,13 lb)

Cable Length 3 m (118,11 in)

Sterilization Autoclavable


Sterilization cycles 100

Compatibility 0350, F790, 00401.10, 00401.20, 5365A, 6429A



pcs: 1

## 00404.10


**Description**

Cable for neutral plate US type

0,06 kg (0,13 lb)

Cable Length 5 m (169,85 in)

Sterilization NO


Sterilization cycles NO

Compatibility 0350, F790, 00401.10, 00401.20, 5365A, 6429A



pcs: 1

## 00404.11


**Description**

Cable for neutral plate US type AUTOCLAVE

0,06 kg (0,13 lb)

Cable Length 5 m (169,85 in)

Sterilization Autoclavable

Sterilization cycles 100


Compatibility 0350, F790, 00401.10, 00401.20, 5365A, 6429A



pcs: 1

<b>Description</b>	Bipolar cable EUR
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	310 Bipolar Forcebs series

00411.00




50 D	80 D
x x + + + +	x + + + +
80	80
# # # #	# # # #
120	120
160	160
200	200
300 HP	TOUCH 200
	FLASH 120
	FLASH 160
	FLASH 200
	300 HP
	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



pcs: 1

<b>Description</b>	Bipolar cable TWIN
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	310 Bipolar Forcebs series

00412.00



50 D	80 D
x x # # # #	x + # # # #
80	80
# # # #	# # # #
120	120
160	160
200	200
300 HP	TOUCH 200
	FLASH 120
	FLASH 160
	FLASH 200
	300 HP
	400 HP
	TOUCH HP
	TOUCH HP (Argon)
	SURTRON EVAC
	SURTRON ABC



pcs: 1

<b>Description</b>	Artery Sealer Bipolar cable
<b>Weight</b>	0,06 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	110-750NS, 110-755NS, 110-760NS,

50 D	80 D
x x # # # #	x + # # # #
80	80
# # # #	# # # #
120	120
160	160
200	200
# TOUCH 200	# FLASH 120
# FLASH 160	# FLASH 200
# 300 HP	# 400 HP
# TOUCH HP	# TOUCH HP (Argon)
# SURTRON EVAC	# SURTRON ABC



pcs: 1


<b>Description</b>	Bipolar cable US 3m
<b>Weight</b>	0,60 kg (0,13 lb)
<b>Cable Length</b>	3 m (118,11 in)
<b>Sterilization</b>	NO
<b>Sterilization cycles</b>	NO
<b>Compatibility</b>	110-700

50 D	80 D
x x # # # #	x + # # # #
80	80
# # # #	# # # #
120	120
160	160
200	200
# TOUCH 200	# FLASH 120
# FLASH 160	# FLASH 200
# 300 HP	# 400 HP
# TOUCH HP	# TOUCH HP (Argon)
# SURTRON EVAC	# SURTRON ABC




pcs: 1

00414.00



## 00415.00

**Description**

Bipolar cable ENDO (MP3 - F3)

0,06 kg (0,13 lb)

3 m (118,11 in)

Autoclavable

100


**Weight****Cable Length****Sterilization****Sterilization cycles**

x	50 D	x	80 D	x	80	x	120	x	160	x	200	x	200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	



pcs: 1

## 00416.00

**Description**

Bipolar cable ENDO (MP3 - F3)

0,06 kg (0,13 lb)

3 m (118,11 in)

Autoclavable

100


**Weight****Cable Length****Sterilization****Sterilization cycles**

x	50 D	x	80 D	x	80	x	120	x	160	x	200	x	200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		



pcs: 1

## 00417.00

**Description**

Bipolar cable ENDO3 (MP2-F4)

0,06 kg (0,13 lb)

3 m (118,11 in)

Autoclavable

100

**Weight****Cable Length****Sterilization****Sterilization cycles**

x	50 D	x	80 D	x	80	x	120	x	160	x	200	x	200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		



pcs: 1

## 00418.00

**Description**

Bipolar cable ENDO3 (MP2-F4)

0,06 kg (0,13 lb)

3 m (118,11 in)

Autoclavable

100

**Weight****Cable Length****Sterilization****Sterilization cycles**

x	50 D	x	80 D	x	80	x	120	x	160	x	200	x	200	x	FLASH 120	x	FLASH 160	x	FLASH 200	x	FLASH 200	x	300 HP	x	400 HP	x	TOUCH HP	x	TOUCH HP (Argon)
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		




pcs: 1




ELEC  
TRODE



00500.00



500500.L11



## 152 Short



LED<sup>®</sup>  
30

**Description** Kit of assorted electrode length 5 cm (10pcs)

## SINGLE ITEM

- 500500.L1: Straight thin wire electrode (5pcs)
  - 500500.L3: Loop electrode ø4 mm (5pcs)
  - 500500.L4: Loop electrode ø6 mm (5pcs)
  - 500500.L7: Drop electrode length (5pcs)
  - 500500.L8: Noose electrode length (5pcs)
  - 500500.L9: Straight ball electrode ø3 mm (5pcs)
  - 500500.L2: Bent thin wire electrode (5pcs)
  - 500500.L5: Drop electrode length (5pcs)
  - 500500.L6: Bent thick wire electrode (5pcs)
  - 500500.L10: Bent5 ball electrode ø3 mm (5pcs)

**Compatibility** 00205.00, 00206.00, 755VL

kit 00500.00: assorted pcs 10  
for each items from L1 to L10: pcs 5

<b>Description</b>	Needles for micro-surgery length 1cm
<b>Weight</b>	0,002 kg (0,004 lb)
<b>Size</b>	12x0,05x0,8 mm (0,47x0,002x0,03 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	NO



pcs: 10

## Description

- 152-110 Blade Electrode 7 cm

+ 50 D	ω 200	ω TOUCH 200	ω FLASH 120	ω FLASH 160	ω FLASH 200	ω 300 HP	ω 400 HP	# TOUCH HP
+ 80 D	+ 80							
+ 120	+ 120							
+ 160	+ 160							
152-112	Blade Curved Electrode 7 cm							
152-120	Needle Curved Electrode 7 cm							
152-122	Needle Electrode 13 cm							
152-130	Ball Electrode ø2 mm 6 cm							
152-132	Ball Curved Electrode ø2 mm 5 cm							
152-140	Ball Electrode ø3 mm 6cm							
152-142	Ball Curved Electrode ø3 mm 6 cm							

152-152 Ball Curved Electrode ø4 mm 6 cm

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	#	TOUCH HP
x	+	+	+	+	+	+	+	+	+	+	+		
152-160													Ball Electrode ø5 mm 6 cm
152-162													Ball Curved Electrode ø5 mm 6cm

**Compatibility** 00205.00, 00206.00, 755VL



**Description** Kit of assorted electrode length 10 cm (10pcs)

**SINGLE ITEM**

- 500500.L/L1: Straight thin wire electrode length (5pcs)
- 500500.L/L3: Loop electrode ø4 mm length (5pcs)
- 500500.L/L4: Loop electrode ø6 mm (5pcs)
- 500500.L/L7: Drop electrode length (5pcs)
- 500500.L/L8: Noose electrode length (5pcs)
- 500500.L/L9: Straight ball electrode ø3 mm (5pcs)
- 500500.L/L2: Bent thin wire electrode (5pcs)
- 500500.L/L5: Bent hook electrode (5pcs)
- 500500.L/L6: Bent thick wire electrode (5pcs)
- 500500.L/L10: Straight ball electrode ø3 mm (5pcs)

**Compatibility** 00205.00, 00206.00, 755VL

00500.00/L



+ 50 D	+ 80 D	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+ + 200	+ + 80										
+ + 120	+ + 160										
+ + 160	+ + 200										

pcs:

kit 00500.00\_L/L: assorted pcs 10  
for each items from L-L1 to L-L10 pcs 5

**Description**



- 152-115 Blade Electrode 16 cm
- 152-125 Needle Electrode 13 cm
- 152-145 Ball Electrode ø3mm 13 cm
- 152-165 Ball Electrode ø5mm 14 cm
- 152-175-10 Loop Electrode 10x10 l.14 cm
- 152-190-13 Loop Electrode 20x13 l.14 cm
- 152-190-20 Loop Electrode 20x20 l.15 cm
- 152-195 Conization Electrode 13 cm

**Compatibility** 00205.00, 00206.00, 755VL

+ 50 D	+ 80 D	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+ + 200	+ + 80										
+ + 120	+ + 160										
+ + 160	+ + 200										

pcs: 5

152 Long





# ADAP TERS



# BI PO LAR

# MO NO PO LAR



<b>Description</b>	US Bi-Turbinate Electrode
<b>Size</b>	18,5x2x2 cm (7,3x0,78x78 in)
<b>Weight</b>	0,07 kg (0,154 in)
<b>Tip Thickness</b>	4,5 mm (1,77 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Compatibility</b>	00414.00

110-700



175 mm



50	80	80	120	200	TOU	FLA	FLA	FLA	300	400	TOU	TOU	SUR	SUR	
x	x	x	x	x	x	x	x	x	x	x	#	#	#	#	pcs: 1

110-750NS



27 cm



50 D	80	80	120	160	200	TOUCH-	FLASH-	FLASH-	FLASH	300 H	400 H	TOUCH	TOUCH	SURF	SURF	TR
x	x	x	x	x	x	x	x	x	x	x	x	#	#	#	#	#

---

<b>Description</b>	Bipolar Forceps Artery Sealer
<b>Size</b>	25,5 cm (10 in)
<b>Size</b>	27 cm (10,62 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Tip thickness</b>	3mm
<b>Compatibility</b>	00413.00

110-755NS




25.5 cm



50 D	80 D	80	120	160	200	TOU	FLA:	FLA:	FLA:	300	400	TOU	TOU	SUR	SUR	
-----																
<b>Description</b>	Bipolar Forceps Artery Sealer															
<b>Size</b>	17 cm (6,7in)															
<b>Size</b>	27 cm (10,62 in)															
<b>Sterilization</b>	Autoclavable															
<b>Sterilization cycles</b>	30															
<b>Tip thickness</b>	2 mm															
<b>Compatibility</b>	00413.00															

110-760NS



17 cm



## 310-110-05



11,5 cm



Tip Shape

**Description**

Bipolar Forceps 11,5cm

0,03 kg (0,6lb)

11,5 cm (4,52 in)

0,5 mm (0,019in)

Autoclavable

Sterilization cycles 30

Compatibility 00411.00, 00412.00



pcs: 1

## 310-112-05



11,5 mm



Tip Shape

**Description**

Curved Bipolar Forceps

0,03 kg (0,6lb)

11,5 cm (4,52 in)

0,5 mm (0,019in)

Autoclavable

Sterilization cycles 30

Compatibility 00411.00, 00412.00



pcs: 1

## 310-140-10



20 cm



Tip Shape

**Description**

Bipolar Forceps

0,03 kg (0,6lb)

20 cm (7,87lb)

1 mm (0,039in)

Autoclavable

Sterilization cycles 30

Compatibility 00411.00, 00412.00



pcs: 1

## 310-140-20



20 cm



Tip Shape

**Description**

Bipolar Forceps

0,03 kg (0,6lb)

20 cm (7,87 lb)

2 mm (0,078 in)

Autoclavable

Sterilization cycles 30

Compatibility 00411.00, 00412.00



pcs: 1

<b>Description</b>	Curved Bipolar Forceps
<b>Weight</b>	0,03Kg (0,6lb)
<b>Length</b>	20 cm (7,87 in)
<b>Tip thickness</b>	1 mm (0,039 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00

310-142-10

50 D  
x x 80 D  
x x + 80  
+ + 120  
+ + 160  
+ + 200  
TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
300 HP  
400 HP  
# # # SURTRON EVAC  
# # # SURTRON ABC



pcs: 1

<b>Description</b>	Curved Bipolar Forceps
<b>Weight</b>	0,03 Kg (0,6 lb)
<b>Length</b>	20 cm (7,87 in)
<b>Tip thickness</b>	2 mm (0,078 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00

310-142-20

50 D  
x x 80 D  
x x + 80  
+ + 120  
+ + 160  
+ + 200  
TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
300 HP  
400 HP  
# # # TOUCH HP (Argon)  
# # # SURTRON EVAC  
# # # SURTRON ABC



pcs: 1

<b>Description</b>	Bipolar Forceps Artery Sealer
<b>Weight</b>	0,03 kg (0,6 lb)
<b>Length</b>	25 cm (10 in)
<b>Tip thickness</b>	3 mm (0,11 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00

310-180-10

50 D  
x x 80 D  
x x + 80  
+ + 120  
+ + 160  
+ + 200  
TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
300 HP  
400 HP  
# # TOUCH HP (Argon)  
# # # SURTRON EVAC  
# # # SURTRON ABC



pcs: 1



<b>Description</b>	Angled Bipolar Forceps 20cm
<b>Weight</b>	0,03 kg (0,6 lb)
<b>Length</b>	17 cm (0,07 in)
<b>Tip thickness</b>	2 mm (0,07 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilization cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00

310-180-20

x 50 D  
x 80 D  
x + 80  
+ + 120  
+ + 160  
+ + 200  
TOUCH 200  
FLASH 120  
FLASH 160  
FLASH 200  
300 HP  
400 HP  
# # TOUCH HP (Argon)  
# # # SURTRON EVAC  
# # # SURTRON ABC



pcs: 1




## 310-182-10



<b>Description</b>	Angled Curved Bipolar Forceps
<b>Weight</b>	0,03 kg (0,6 lb)
<b>Length</b>	20 cm (7,87 lb)
<b>Tip thickness</b>	1 mm (0,039 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00



## 310-185-10



<b>Description</b>	Angled Curved Bipolar Forceps
<b>Weight</b>	0,03 kg (0,6lb)
<b>Length</b>	20 cm (7,87 lb)
<b>Tip thickness</b>	1 mm (0,039 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00



## 310-510



<b>Description</b>	Bipolar Electrode
<b>Weight</b>	0,03 kg (0,6 lb)
<b>Length</b>	20 cm (7,87 lb)
<b>Tip thickness</b>	1 mm (0,039 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00



## 310-550



<b>Description</b>	Bipolar Electrode
<b>Weight</b>	0,03 kg (0,06lb)
<b>Length</b>	20 cm (7,87lb)
<b>Tip thickness</b>	0,5 mm (0,02in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00



310-590

<b>Description</b>	Bipolar Forcep 20 cm - curved 2
<b>Weight</b>	0,03 Kg (0,6 lb)
<b>Length</b>	18 cm (7 lb)
<b>Tip thickness</b>	0,5 mm (0,02 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00411.00, 00412.00

50 D	80 D	120	160	200
x x + + + +	+ + + + + +	# # # # # #		
TOUCH 200	FLASH 120	FLASH 160	FLASH 200	



pcs: 1



18 cm

Tip Shape



<b>Description</b>	Bipolar Forcep Clamp Scissors
<b>Weight</b>	0,05 kg (0,1 lb)
<b>Length</b>	15 cm (5,9 lb)
<b>Tip thickness</b>	0,5 mm (0,02 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00418.00

50 D	80 D	120	160	200
x x x x x x	+ + + + + +	# # # # # #		
TOUCH 200	FLASH 120	FLASH 160	FLASH 200	



pcs: 1



15 cm

410-100-15

<b>Description</b>	Bipolar Forcep Clamp Scissors
<b>Weight</b>	0,05 kg (0,1 lb)
<b>Length</b>	19 cm (7,5 lb)
<b>Tip thickness</b>	0,5 mm (0,02 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00418.00

50 D	80 D	120	160	200
x x x x x x	+ + + + + +	# # # # # #		
TOUCH 200	FLASH 120	FLASH 160	FLASH 200	



pcs: 1



19 cm

410-100-19

<b>Description</b>	Bipolar Forcep Scissors
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Length</b>	18 cm (7,0 lb)
<b>Tip thickness</b>	2 mm (0,078 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00418.00

50 D	80 D	120	160	200
x x + + + +	+ + + + + +	# # # # # #		
TOUCH 200	FLASH 120	FLASH 160	FLASH 200	



pcs: 1




18 cm

Tip Shape



## 410-200-21



21 cm



Tip Shape


<b>Description</b>	Bipolar Forcep Scissors
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Length</b>	21 cm (7,0 lb)
<b>Tip thickness</b>	2 mm (0,078 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00418.00

x	50 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80	+	+	+	#	#	SURTRON EVAC
x	120	+	+	+	#	#	SURTRON ABC
x	160	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC



pcs: 1

## 410-200-23



23 cm



Tip Shape


<b>Description</b>	Bipolar Forcep Scissors
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Length</b>	23 cm (7,0 lb)
<b>Tip thickness</b>	2 mm (0,078 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	00418.00

x	50 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80	+	+	+	#	#	SURTRON EVAC
x	120	+	+	+	#	#	SURTRON ABC
x	160	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC



pcs: 1

## 330-160



18 cm

Tip Shape

<b>Description</b>	Monopolar Forceps 21cm
<b>Weight</b>	0,05Kg (0,11lb)
<b>Length</b>	18 cm (7,0 lb)
<b>Cutter length</b>	45 mm (1,7 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	0040202

x	50 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80	+	+	+	#	#	SURTRON EVAC
x	120	+	+	+	#	#	SURTRON ABC
x	160	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC



pcs: 1

## 330-134-20



21 cm

Tip Shape

<b>Description</b>	Monopolar Surgical Scissors
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Length</b>	21 cm (8,2 lb)
<b>Tip thickness</b>	2 mm (0,078 in)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	0040202

x	50 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80 D	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP (Argon)
x	80	+	+	+	#	#	SURTRON EVAC
x	120	+	+	+	#	#	SURTRON ABC
x	160	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC
x	200	+	+	+	#	#	SURTRON ABC



pcs: 1

# EVA CUA TION KIT



## 00900.01/12



STERILE EO

<b>Description</b>	Handle evacuation kit for surgical smoke evacuator
<b>Pipe Diameter</b>	ø22 cm (3,93 in)
<b>Pipe Length</b>	18,5 cm (7,2 in)
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Sterilization</b>	NO
<b>Sterilisation cycles</b>	NO

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x	x	x	x	x	x	x	x	x	x	x	x	x	x	+	x

pcs: 12

## 00900.02/06



<b>Description</b>	Evacuation kit for surgical smoke evacuator
<b>Pipe Diameter</b>	ø22 cm (0,86 in)
<b>Pipe Length</b>	240 cm (94,48 in)
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Sterilization</b>	NO
<b>Sterilisation cycles</b>	NO
<b>Compatibility</b>	00900.01/12, 00900.05, 00900.07, 00900.08/12,

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x	x	x	x	x	x	x	x	x	x	x	x	x	x	+	x

pcs: 6

## 00900.04



134°C  
 CICLI 100 CYCLES

<b>Description</b>	Autoclavable evacuation tube D.22
<b>Pipe Diameter</b>	ø7 mm (0,27 in)
<b>Pipe Length</b>	2,40 m (94,5 in)
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	100
<b>Compatibility</b>	00900.01/12, 00900.05, 00900.06/06, 00900.07, 00900.08/12

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
x	x	x	x	x	x	x	x	x	x	x	x	x	x	+	x

pcs: 1

# 00900.05

<b>Description</b>	Autoclavable evacuation tube D.22
<b>Pipe Diameter</b>	ø7 mm (0,27 in)
<b>Pipe Length</b>	2,40 m (94,5 in)
<b>Weight</b>	0,05 kg (0,11lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	100
<b>Compatibility</b>	00900.02, 00900.04

50 D	FLASH 120
x x 80 D	FLASH 160
80	FLASH 200
x x x x x x	TOUCH 200
120	TOUCH HP
160	TOUCH HP Argon
200	SURTRON EVAC
	SURTRON ABC



pcs: 1



ø7 mm

ADAPTER  
7/22

<b>Description</b>	Bacterial/Viral Filter D.22
<b>Pipe Diameter</b>	ø22 mm (0,867 in)
<b>Connection Diameter</b>	ø15 mm (0,59 in)
<b>Sterilization</b>	NO
<b>Sterilisation cycles</b>	NO
<b>Compatibility</b>	00900.04, 00900.02/06, 00900.FU01/10

50 D	FLASH 120
x x 80 D	FLASH 160
80	FLASH 200
x x x x x x	TOUCH 200
120	TOUCH HP
160	TOUCH HP (Argon)
200	SURTRON EVAC
	SURTRON ABC



pcs: 1



ø22 mm

<b>Description</b>	Autoclavable handle evacuation kit
<b>Pipe Diameter</b>	ø22 cm (0,86 in)
<b>Pipe Length</b>	150 cm (59 in)
<b>Weight</b>	0,05 kg (0,11lb)
<b>Sterilization</b>	Autoclavable
<b>Sterilisation cycles</b>	30
<b>Compatibility</b>	755VL, 00205.00/40, 00206.00/40 00900.02/06, 00900.04

50 D	FLASH 120
x x 80 D	FLASH 160
80	FLASH 200
x x x x x x	TOUCH 200
120	TOUCH HP
160	TOUCH HP (Argon)
200	SURTRON EVAC
	SURTRON ABC



pcs: 1

# 00900.07



## 00900.08/12



**STERILE EO**

<b>Description</b>	Evacuation kit for TROCAR Luer Lock
<b>Pipe Diameter</b>	ø22 mm (0,86 in)
<b>Pipe Length</b>	240 cm (94,48 in)
<b>Weight</b>	0,05 kg (0,11 lb)
<b>Sterilization</b>	NO
<b>Sterilisation cycles</b>	NO
<b>Compatibility</b>	00900.02, 00900.04

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	
FLASH 120		FLASH 160		FLASH 200		300 HP		400 HP		TOUCH HP		TOUCH HP (Argon)	
SURTRON EVAC		SURTRON ABC											



pcs: 12

## 00900.FU01/10\*



<b>Description</b>	ULPA Filter for SURTRON EVAC - red
<b>Filter Diameter</b>	22 mm 7/8" (0,867 in)
<b>Particles diameter</b>	0,12 micron
<b>Filter Type</b>	ULPA with active carbon
<b>Filter efficiency</b>	99,999X %
<b>Filter duration</b>	00900.FST/10: 5 hours
<b>Filter duration</b>	00900.FU01/10: 20 hours

50 D	x	80 D	x	80	x	120	x	160	x	200	x	TOUCH 200	
FLASH 120		FLASH 160		FLASH 200		300 HP		400 HP		TOUCH HP		TOUCH HP (Argon)	
SURTRON EVAC		SURTRON ABC											




pcs: 3

# TROL LEY



## TR003



**Description**  
**Size WxHxD**  
**Weiht**  
**Shelf size**  
**Color**  
**Frame Material**

**Wheel type**

Trolley 3 shelves

80x46x46 cm (31x18,1x18,1 in)

18 kg (34,47 lb)

43x36 cm (16,9x14,17 in)

RAL7035 (Light Grey)

Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.

80 mm pirouette, anterior wheels whit block system.

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	x x	x x	x x	x x	+ +	

pcs: 1

WxHxD: 13x47x67 cm

Weight: 18 kg

## TR003W



**Description**  
**Size WxHxD**  
**Weiht**  
**shelf size**  
**Color**  
**Frame Material**

**Wheel type**

Trolley 3 shelves

80x54x46 cm (31x21,2x18,1 in)

18,5 kg (34,47 lb)

40x40 cm (15,7x15,74 in)

RAL7035 (Light Grey)

Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.

80 mm pirouette, anterior wheels, whit block system.

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	x x	x x	x x	x x	+ +	



pcs: 1

WxHxD: 13x47x67 cm

Weight: 18,5 kg

## TR004



**Description**  
**Size WxHxD**  
**Weiht**  
**Shelf size**  
**Color**  
**Frame Material**

**Wheel type**

Trolley 4 shelves

90x46x46 cm (35,43x18,1x18,1 in)

20 kg (39,6 lb)

40x36 cm (15,7x14,17 in)

RAL7035 (Light Grey)

Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.

80 mm pirouette, anterior wheels whit block system.

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP	TOUCH HP (Argon)	SURTRON EVAC	SURTRON ABC
+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	+ +	x x	x x	x x	x x	+ +	



pcs: 1

WxHxD: 13x47x75 cm

Weight: 20 kg

# TR005

<b>Description</b>	Trolley 5 shelves
<b>Size WxHxD</b>	112x46x46 cm (44x18,1x18,1 in)
<b>Weight</b>	19 kg (41,9 lb)
<b>Shelf size</b>	43x36 cm (16,9x14,17 in)
<b>Color</b>	RAL7035 (Light Grey)
<b>Frame Material</b>	Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
<b>Wheel type</b>	80 mm pirouette, anterior wheels with block system.

50 D	FLASH 120
+ 80 D	FLASH 160
+ 80	FLASH 200
+ 120	
+ 160	
+ 200	
TOUCH 200	

pcs: 1  
► shelves WxHxD: 13x47x67 cm  
 side columns WxHxD: 104x15x20 cm  
 Weight: 27,5 kg



<b>Description</b>	Trolley 5 shelves
<b>Size WxHxD</b>	112x46x54 cm (44x18,1x21,2 in)
<b>Weight</b>	24,5 kg (54 lb)
<b>Shelf size</b>	43x36 cm (16,9x14,17 in)
<b>Color</b>	RAL7035 (Light Grey)
<b>Frame Material</b>	Trolley made of aluminum with adjustable shelves made with steel, covered with a coat of paint.
<b>Wheel type</b>	80 mm pirouette, anterior wheels with block system.

50 D	FLASH 120
+ 80 D	FLASH 160
+ 80	FLASH 200
+ 120	
+ 160	
+ 200	
TOUCH 200	

pcs: 1  
► shelves WxHxD: 13x47x67 cm  
 side columns WxHxD: 104x15x20 cm  
 Weight: 27,5 kg

# TR005W




<b>Description</b>	Support for Cylinder
<b>Size WxHxD</b>	18x6x20 cm (7x2,3x7,87 in)
<b>Weight</b>	0,5 kg (1,1 lb)
<b>Color</b>	RAL7016 (Anthracite Gray)
<b>Frame Material</b>	Shelves made with steel, covered with a coat of paint.

50 D	FLASH 120
+ 80 D	FLASH 160
+ 80	FLASH 200
+ 120	
+ 160	
+ 200	
TOUCH 200	

pcs: 1  
 WxHxD: 67x47x23 cm

# TRCYLSUP01



## TRDRAWERO1



**Description**  
**Size WxHxD**

Drawer for trolley type TR W  
38,5x10x28,7 cm (15x3,9x11,3 in)

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP Argon	SURTRON EVAC	SURTRON ABC
+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	+ + +	x x x	x x x	x x x



pcs: 1

## TRDRAWERO1W



**Description**  
**Size WxHxD**


Drawer for trolley type TR W  
38,5x10x28,7 cm (15x3,9x11,3 in)

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP Argon	SURTRON EVAC	SURTRON ABC
x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x



pcs: 1

## TRSTO01



**Description**  
**Size WxHxD**

**Weight**  
**shelf size**  
**Color**

**Frame Material**

**Wheel type**

SURTRON TOUCH HP Trolley  
53x108,5x57,5 cm (20x42,7x22,6 in)  
33 Kg (72,75lb)  
43,5x10x38 mm (17x3,9x15 in)

RAL7016 (Anthracite gray)  
RAL6027 (Light Green)  
RAL9005 (Intense Black)

Trolley with column made of aluminum with adjustable shelves made with steel, covered with a coat of paint.

80 mm pirouette, anterior wheels with block system.

50 D	80 D	80	120	160	200	TOUCH 200	FLASH 120	FLASH 160	FLASH 200	300 HP	400 HP	TOUCH HP Argon	SURTRON EVAC	SURTRON ABC
x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x	x x x



pcs: 1

WxHxD: 100x60x18 cm    Weight: 33 kg

*TRSTBO1*

<b>Description</b>	Basket for ST Trolley
<b>Size WxHxD</b>	53x108,5x57,5 cm (20x42,7x22,6 in)
<b>Weiht</b>	3,2 kg (7 lb)
<b>shelf size</b>	43,5x10x38 mm (17x3,9x15 in)
<b>Color</b>	RAL9005 (Intense Black)
<b>Frame Material</b>	Basket made with steel, covered with a coat of paint.



50 D	TOUCH 200
80 D	FLASH 120
80	FLASH 160
120	FLASH 200
160	300 HP
200	400 HP
	TOUCH HP
	TOUCH HP Argon
	SURTRON EVAC
	SURTRON ABC



pcs: 1

<b>Description</b>	Support section for ST Trolley
<b>Size WxHxD</b>	43,5x13,5x38 cm (17x5,3x15 in)
<b>Weiht</b>	5,75Kg (2,26lb)
<b>Color</b>	RAL6027 (Light Green)
<b>Frame Material</b>	Holder made with steel, covered with a coat of paint.

50 D	TOUCH 200
80 D	FLASH 120
80	FLASH 160
120	FLASH 200
160	300 HP
200	400 HP
	TOUCH HP
	TOUCH HP Argon
	SURTRON EVAC
	SURTRON ABC



pcs: 1

<b>Description</b>	Basket for ST Trolley
<b>Size WxHxD</b>	14,5x30x19,3 cm (5,7x11,8x7,6 in)
<b>Weiht</b>	2,4 Kg (5,29 lb)
<b>Color</b>	RAL6027 (Light Green)
<b>Frame Material</b>	Basket holder applicable to the Surtron Touch HP Trolley made with steel, covered with a coat of paint.

50 D	TOUCH 200
80 D	FLASH 120
80	FLASH 160
120	FLASH 200
160	300 HP
200	400 HP
	TOUCH HP
	TOUCH HP Argon
	SURTRON EVAC
	SURTRON ABC



pcs: 1

*TRSTCS1*



## *NOTE*





## LED SpA

Via Selciatella, 40  
04011 - Aprilia (LT) Italy

[www.led.it](http://www.led.it) - [info@led.it](mailto:info@led.it)

Tel. +39 06 92870045  
Fax. +39 06 92870046

