

# COVER-ENT®

## STERILE SINGLE-USE SHEATHS FOR FLEXIBLE NASOPHARYNGOSCOPES



Made  
in Italy



Quality Control  
on 100% of  
production



+ Maintains the quality of  
vision

Sterile, practical, and  
ready to use

Greater safety and  
productivity

Effective barrier against  
microorganisms

Lower risks of  
infection and  
adverse events

Preserving the efficiency  
of nasopharyngoscopes

Reducing maintenance  
costs and times

# COVER

CE 0477

Regulation (EU) 2017/745

Cover S.r.l, founded in 2005, was born from the experience of entrepreneurs working in specific sectors since the 80's with the mission to develop and produce high quality medical devices. The goal has been achieved thanks to our long experience and following our clients' suggestions.

# COVER-ENT®

## STERILE SINGLE-USE SHEATHS FOR FLEXIBLE NASOPHARYNGOSCOPES



### PACKAGING:

20 pcs per package,  
5 packs per box

WE ARE ABLE TO PRODUCE SHEATHS FOR ANY  
TYPE OF NASOPHARYNGOSCOPE NOT IN THE  
BROCHURE



Certified Quality Management System  
UNI CEI EN ISO  
13485:2021



Distributed by:

# ENDOVISION

SURGICAL INSTRUMENTS / EQUIPMENT / CONSUMABLES

Manufactured by:

# COVER

A YUSHU Healthcare Company

YUSHU  
HEALTHCARE



ENDOSCOPE BRAND	ENDOSCOPE MODEL	COVER SHEATH CODE
OLYMPUS	ENF-XP	GNF20X30
	ENF-P4	GNF42X300
	ENF-GP	GNF42X300
	ENF-GP2	GNF33X300
	ENF-V	GNF44X37
	ENF-V2	GNF42X300
	ENF-V3	GNF30X30P
	ENF-V4	GNF30X32O
	ENF-VH	GNF4230L3
	ENF-VH2	GNF4230L3
	ENF-VQ	GNF4230L3
STORZ	11101SK-SK2	GNF30X27
	11101RP1-2	GNF39X30S
	11101CM	GNF39X30
	11101HD	GNF39X30/39X30V
	11101VP	GNF39X30V
	11102CM	GNF30X30/30X30L
	11102CM	GNF30X30/30X30L
PENTAX	FNL-7RP3	GNF30X30S
	FNL-10RP3	GNF42X300
	FNL-10RBS	GNF42X300
	VNL8-J10	GNF30X30L
	VNL11-J10	GNF4230L3
	VNL1170K	GNF4430L2
	VNL1190K	GNF4430L2
XION	EF-N SLIM	GNF30X32
	EF-N / EF-NS	GNF39X32
	EV-N	GNF44X32
	EV-NC	GNL39X30S
	EV-NE	GNF4430X2
	XN-HD	GNF4232L3
	XN-HD	GNF4232L3
FENTEX	700710FX	GNF35X30
MACHIDA	ENT-30PIII	GNF35X30M
OTOPRONT	VN-S	GNF39X28
ORL VISION	RS1-RX1	GNFORLV01
OPTIM	Entity	GNF39X30E
HENKE SASS WOLF	8300019728	GNF35X30

ENT 07-25

