

# COVER-ENT<sup>®</sup> SHEATHS FOR RIGID ENDOSCOPES

*According to Prof. Paolo Castelnovo*



Made  
in Italy



Quality Control  
on 100% of  
production



+ It does not increase patient  
discomfort during the endoscopic  
examination

Maintains endoscopic image  
quality

It allows you to preserve the  
sterility of the telescope by being  
able to visit more patients without  
waiting times

Extends the life of the telescope by  
reducing sterilization cycles

Effective barrier against  
microorganism

COVER

CE 0477

Regulation (EU) 2017/745

Cover S.r.l, founded in 2005, was born from the experience of entrepreneurs working in specific sectors since the 80's with the mission to develop and produce high quality medical devices. The goal has been achieved thanks to our long experience and following our clients' suggestions.

# STERILE SINGLE-USE COVER-ENT<sup>®</sup> SHEATHS FOR RIGID ENDOSCOPES



## PACKAGING:

20 pcs per package,  
5 packs per box



FOR RIGID ENDOSCOPES WITH CONNECTION TYPE STORZ,  
COLLIN HENKE SASS WOLF AND XION

CODE	DIRECTION OF VIEW	OPTIC DIAMETER IN MM	WORKING LENGTH IN MM
GFR27110	0°	2,7	110
GFR271130	30°	2,7	110
GFR30110	0°	3	110
GFR30140	0°	3	140-150
GFR301430	30°	3	140-150
GFR301445	45°	3	140-150
GFR27180	0°	2,7	180
GFR271830	30°	2,7	180
GFR271845	45°	2,7	180
GFR40180	0°	4	170-180
GFR401830	30°	4	170-180
GFR401845	45°	4	170-180
GFR401870	70°	4	170-180

ENT 07-24

Certified Quality Management System  
UNI CEI EN ISO  
13485:2021



Distributed by:

**ENDOVISION**  
SURGICAL INSTRUMENTS / EQUIPMENT / CONSUMABLES

Manufactured by:

**COVER**

A YUSHU Healthcare Company

**YUSHU**  
HEALTHCARE



C.so Re Umberto, 8 , 10121 Torino (Italy)  
Tel **+39 011 0438523**  
email: [info@yushuhealthcare.com](mailto:info@yushuhealthcare.com)  
[www.yushuhealthcare.com](http://www.yushuhealthcare.com)

